Minuteman® G1 Interspinous Interlaminar Fusion Device
Minuteman® HA-Coated G1 Interspinous Interlaminar Fusion Device

IMPORTANT INFORMATION ON THE MINUTEMAN INTERSPINOUS INTER-LAMINAR FUSION DEVICE
WICHTIGE INFORMATIONEN ZUM MINUTEMAN INTERSPINOSEN INTERLAMINAREN FUSIONSIMPLANTAT
INFORMACIÓN IMPORTANTE ACERCA DEL DISPOSITIVO MINUTEMAN DE FUSIÓN INTERLAMINAR E

The surgeon must be fully conversant with all aspects of the surgical technique and know the indicator and contra-indicators of this type of implant. Before beginning the surgical procedure, the surgeon must be acquainted with the specific technique for insertion of the implant, which is available from the manufacturer.

Use only new implants for each case. Implants are not to be re-sterilized. Re-use of the implants could result in infection.

General Conditions of Use:

The information contained in this package insert is necessary but not sufficient for the use of this device. This information is not intended as a substitute for the professional judgment, skill and experience of the surgeon in: careful patient selection; preoperative planning; device selection; knowledge of the anatomy and biomechanics of the spine;

Intraoperative:

- As part of the preoperative examination, the surgeon must check that no factors, especially biological and biomechanical, will affect the correct performance of the implant during the operation and postoperative period.
Only patients that meet the criteria described in the indications For Use should be selected.
Patient conditions and/or pre-dispositions such as those addressed in the aforementioned contraindications should not be selected.

Postoperative:

- At all times, extreme caution should be used around the spinal cord and nerve roots. Damage to nerves may result in loss of neurological functions.
Breakage, slippage or misuse of instruments or implant components may cause injury to the patient or operative personnel.
The correct selection of the type and size of implant appropriate to the patient, and positioning of the implant, are extremely important to allow maximum chances for a successful surgical result.

Contraindications:

- The physician's postoperative directions and warnings to the corresponding patient are important to allow maximum chances for a successful surgical result. Detailed instructions on the use and limitations of the device should be given to the patient.
If a nonunion develops or if the components loosen, bend, and/or break, the device(s) should be revised and/or removed immediately before serious injury occurs.

Sterility:

Gamma irradiation is indicated by the "Sterile-R" symbol on the Minuteman implant labeling. These devices remain sterile as long as the package integrity has not been violated.

Packaging:

Packages for each of the Minuteman components should be intact upon receipt. Inspect each package prior to use and do not use the component if any seal or cavity is damaged or breached or if the expiration date has been exceeded.

Handling and storage:

The Minuteman implants and instruments must be stored with care. Before use, inspect all instrumentation for proper function, possible damage, wear or non-function. Damaged or defective instruments should not be used.

Recommendations for Care and Cleaning of Minuteman Surgical Instrument System:

Both physical and chemical (detergent) processes are necessary to minimize the bioburden on all soiled items. Chemical (detergent) cleaners alone cannot remove all soil and debris, therefore a careful manual cleaning of each item is essential for maximum decontamination.

Sonication:

- Remove the instruments from bath and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumen, holes, and difficult-to-reach areas.
Place prepared cleaning agents in a sonication unit. Completely submerge instruments in cleaning solution and sonicate for 10 minutes at 45-50 kHz.

Automated Washer:

- Remove instruments from bath and rinse in purified water for a minimum of 1 minute. Thoroughly and aggressively flush lumen, holes, and difficult-to-reach areas.
Place instruments in a suitable washer/disinfecter basket and process through a standard washer/disinfecter cleaning cycle.

USA Contact

Spinal Simplicity, LLC
10995 Quivira Road
Overland Park, Kansas
66210 USA
T: 913.451.4414
information@spinalsimplicity.com

Europe Contact
Maastrichtlerlaan 127-129
6291 EN Vaais
T: +31 (0)43 306 33 20
F: +31 (0)43 306 33 38
info@meditech.nl
K.v.K nr.: 1408.6283

ENGLISH

IMPORTANT INFORMATION ON THE MINUTEMAN INTERSPINOUS INTER-LAMINAR FUSION DEVICE

Purpose:

The Minuteman System is an interspinous interlaminar fusion device intended for the stabilization of the non-cervical spine (T1-S1). When using the percutaneous technique at levels of L3 or above, a preoperative CT is recommended to ensure a safe trajectory to the spine.

Description:

The Minuteman consists of spinous process plates and surgical instruments. Implants are made of Titanium alloy Ti-6AL-4V ELI per ASTM F136 standards and Ti-6AL-4V per ASTM F1472 standards. Caution: The Minuteman is manufactured from a titanium alloy which is known to produce MRI artifacts.

Indication for Use:

The Minuteman System is an interspinous interlaminar fusion device intended to provide stabilization of the non-cervical spine (T1-S1). It can be used as an adjunct to interbody and/or posterior fusion or as a standalone fixation device.

- Degenerative disc disease
Spinal stenosis
Spondylolisthesis

The Minuteman can be used in single or multilevel constructs. It can be used alone or in combination with other spinal devices (e.g. pedicle screws, facet screws, anterior plates or interbody cages). The Minuteman can be used with or without bone grafting.

Contraindications:

- The Minuteman is contraindicated in patients with:
An allergy to titanium or titanium alloy
Spinal anatomy or disease that would prevent implantation of the device, or cause the device to be unstable in situ such as:
Significant instability of the lumbar spine;
An ankylosed segment at the affected level(s);
Acute fracture of the spinous process or pars interarticularis;
Significant scoliosis (Cobb angle greater than 25 degrees);
Diagnosis of severe osteoporosis
Active systemic infection or infection localized to the site of implantation
Insufficient quality or quantity of bone, which would inhibit rigid device fixation
Incompetent or missing posterior arch
Pregnancy

Possible Adverse Effects:

- Implant dislodgement/migration
Implant not positioned correctly
Fracture of the spinous process
Additional surgery, including removal of the Minuteman device
Foreign body reaction
Mechanical failure of the device
Failure of the device/procedure to improve symptoms and/or function

Surgery related:

- Myocardial infarction
Infection
Blood vessel damage/bleeding
Deep vein thrombosis
Hematoma
Pneumonia
Neurological system compromise
Stroke
Nerve injury or spinal cord damage
Paralysis
Thrombus formation
Graft donor site complications, including pain, fracture or wound healing complications
Dural tear
Wound dehiscence or delayed healing
Pain/discomfort at the operative site
Death

Warnings:

The Minuteman components are not to be connected to the components of devices from another manufacturer. Titanium implants should not be mixed with stainless steel implants in the same construct. A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery, where many extenuating circumstances may compromise the results.

Precautions:

Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of implants, are important considerations in the successful utilization of the Minuteman by the surgeon. Further, the proper selection and compliance of the patient will greatly affect the results. Patients who smoke have been shown to have an increased incidence of non-unions. These patients should be advised of this fact and warned of the consequences. Obese, malnourished, and/or alcoholic patients are poor candidates for spinal fusion.

- e. 2 minute detergent wash with hot tap water (64-66°C/146-150°F)
- f. 15 second hot tap water rinse
- g. 2 minute thermal rinse (80-93°C/176-200°F)
- h. 10 second purified water rinse (64-66°C/146-150°F)
- i. 7 to 30 minute hot air dry (116°C/240°F)

NOTE: The washer/disinfectant manufacturer's instructions should be strictly adhered to.

3. Visually inspect to determine if all visible soil has been removed from the surfaces, lumen, cannulae, crevices, serrations, threading, etc. If visible soil remains, repeat the cleaning/disinfecting procedure.

Caution: certain cleaning solutions such as those containing formalin, glutaraldehyde, mercury, active chlorine, chloride, bromine, bromide, iodine, iodide, and/or alkaline cleaners may damage some instruments. Such cleaning solutions should not be used.

#### Recommendations for Sterilization of Minuteman Surgical Instrument System:

For typical steam autoclave cycles, the following are recommended times and temperatures:

##### 1. Prevacuum Sterilizer:

Wrapped cases, trays and instruments, or cases, trays and instruments should be exposed to 135° C (275° F) for at least 3 minutes. Dry for 30 minutes.

#### Product Complaints:

Communicate suspected deficiencies in product quality, identity, durability, reliability, safety, effectiveness and/or performance directly to Spinal Simplicity, LLC. When filing a complaint provide component name(s), part number(s), lot number(s), your name and address, the nature of the complaint and patient case number. Sterilize and return all component(s) to Spinal Simplicity, LLC. Notify Spinal Simplicity, LLC immediately of any complaints and any incidents resulting in patient death or serious injury.

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