Minuteman® G1 Interspinous Interlaminar Fusion Device  
 Minuteman® G3 Interspinous Interlaminar Fusion Device  
 Minuteman® HA-Coated G1 Interspinous Interlaminar Fusion Device

IMPORTANT INFORMATION ON THE MINUTEMAN INTERSPINOUS INTER-LAMINAR FUSION DEVICE  
 WICHTIGE INFORMATIONEN ZUM MINUTEMAN INTERSPINOSEM INTERLAMINAREN FUSIONSIMPLANTAT  
 INFORMACION IMPORTANTE ACERCA DEL DISPOSITIVO MINUTEMAN DE FUSION INTERLAMINAR E  
 INTERSPINOZA INFORMACION IMPORTANTE CONCERNANT LE DISPOSITIF DE FUSION INTERLAMINARE  
 MINUTEMAN INFORMAZIONI IMPORTANTI SUL DISPOSITIVO DI FUSIONE INTERSPINOZA INFORMAZIONI  
 MINUTEMAN TÄRKEÄÄ TIETOA SELKÄRANGAN NIKAMAVÄLILEYNNISÄISESTÄ MINUTEMAN-  
 FUUSILOITTEESTA VIKTIG INFORMATION OM MINUTEMAN-ANORDNING FÖR BAKRE FUSION MELLAN  
 SPINALUTSKOTT/LAMINAE BELANGRIKE INFORMATIE OVER HET MINUTEMAN INTERSPINAAL  
 INTERLAMINAR FUSIEHULPMIDDEL WAZNE INFORMACJE DOTYCZĄCE MIĘDZYWARSTWOWEGO  
 IMPLANTU MINUTEMAN DO STABILIZACJI MIĘDZYWYROSTKOWEJ INFORMAÇÕES IMPORTANTES SOBRE O  
 DISPOSITIVO MINUTEMAN DE FUSÃO INTERLAMINAR PARA APÓFISES ESPINHOSES INFORMAÇÕES  
 IMPORTANTES SOBRE O DISPOSITIVO MINUTEMAN DE FUSÃO INTERLAMINAR PARA PROCESSOS  
 INTERESPINHOSES MINUTEMAN INTERSPINOZ INTERLAMINAR FUZYON CİHAZI HAKKINDA ÖNEMLİ  
 BİLGİLER

**USA Contact**  
 Spinal Simplicity, LLC  
 10995 Quivira Road  
 Overland Park, Kansas  
 66210 USA  
 T: 913.451.4414  
 information@spinalsimplicity.com

**EC REP**  
**Europe Contact**  
 Maastrichterlaan 127-129  
 6291 EN Vaais  
 T: +31 (0)43 306 33 20  
 F: +31 (0)43 306 33 38  
 info@meditech.nl  
 K.v.K nr.: 1408.6283

**ENGLISH** IMPORTANT INFORMATION ON THE MINUTEMAN INTERSPINOUS INTER-LAMINAR FUSION DEVICE

**Purpose:**  
 The Minuteman System is an interspinous interlaminar fusion device intended for the stabilization of the non-cervical spine (T1-S1). When using the percutaneous technique at levels of L3 or above, a preoperative CT is recommended to ensure a safe trajectory to the spine. The percutaneous technique may not be performed for levels above L1.

**Description:**  
 The Minuteman consists of spinous process plates and surgical instruments. Implants are made of Titanium alloy Ti-6AL-4V ELI per ASTM F136 standards and Ti-6AL-4V per ASTM F1472 standards.  
 Caution: The Minuteman is manufactured from a titanium alloy which is known to produce MRI artifacts. Patients should be advised to disclose the presence of the Minuteman prior to an MRI exam. Failure to do so may affect the quality of diagnostic information obtained from these scans. The Minuteman is MRI safe.

**Indication for Use:**  
 The Minuteman System is an interspinous interlaminar fusion device intended to provide stabilization of the non-cervical spine (T1-S1). It can be used as an adjunct to interbody and/or posterior fusion or as a standalone fixation device. The Minuteman is indicated in the treatment of the following conditions:

- Degenerative disc disease
- Spinal stenosis
- Spondyloisthesis

The Minuteman can be used in single or multilevel constructs. It can be used alone or in combination with other spinal devices (e.g. pedicle screws, facet screws, anterior plates or interbody cages). The Minuteman can be used with or without bone grafting.

**Contraindications:**

- The Minuteman is contraindicated in patients with:
- An allergy to titanium or titanium alloy
  - Spinal anatomy or disease that would prevent implantation of the device, or cause the device to be unstable in situ such as
    - Significant instability of the lumbar spine;
    - An ankylosed segment at the affected level(s);
    - Acute fracture of the spinous process or pars interarticularis;
    - Significant scoliosis (Cobb angle greater than 25 degrees);
  - Diagnosis of severe osteoporosis
  - Active systemic infection or infection localized to the site of implantation
  - Insufficient quality or quantity of bone, which would inhibit rigid device fixation
  - Incompetent or missing posterior arch
  - Pregnancy

**Possible Adverse Effects:**

- Device related:**
- Implant dislodgement/migration
  - Implant not positioned correctly
  - Fracture of the spinous process
  - Additional surgery, including removal of the Minuteman device
  - Foreign body reaction
  - Mechanical failure of the device
  - Failure of the device/procedure to improve symptoms and/or function
- Surgery related:**
- Myocardial infarction
  - Infection
  - Blood vessel damage/bleeding
  - Deep vein thrombosis
  - Hematoma
  - Pneumonia
  - Neurological system compromise
  - Stroke
  - Nerve injury or spinal cord damage
  - Paralysis
  - Thrombus formation
  - Graft donor site complications, including pain, fracture or wound healing complications
  - Dural tear
  - Wound dehiscence or delayed healing
  - Pain/discomfort at the operative site
  - Death

**Warnings:**  
 The Minuteman components are not to be connected to the components of devices from another manufacturer. Titanium implants should not be mixed with stainless steel implants in the same construct. A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery, where many extenuating circumstances may compromise the results.

**Precautions:**  
 Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of implants, are important considerations in the successful utilization of the Minuteman by the surgeon. Further, the proper selection and compliance of the patient will greatly affect the results. Patients who smoke have been shown to have an increased incidence of non-unions. These patients should be advised of this fact and warned of the consequences. Obese, malnourished, and/or alcoholic patients are poor candidates for spinal fusion. Patients

with poor muscle and bone quality and/or nerve paralysis are also poor candidates for spinal fusion.

The surgeon must be fully conversant with all aspects of the surgical technique and know the indicator and contra-indicators of this type of implant. Before beginning the surgical procedure, the surgeon must be acquainted with the specific technique for insertion of the implant, which is available from the manufacturer. Never re-use any implant even if it appears unmarked or undamaged. Any implant implanted and then removed must be discarded.

Use only new implants for each case. Implants are not to be re-sterilized. Re-use of the implants could result in infection.

**General Conditions of Use:**

The information contained in this package insert is necessary but not sufficient for the use of this device. This information is not intended as a substitute for the professional judgment, skill and experience of the surgeon in: careful patient selection; preoperative planning; device selection; knowledge of the anatomy and biomechanics of the spine; understanding of the material and the mechanical characteristics of the implants used; training and skill in both spinal surgery and use of associated instruments for implantation; securing the patients cooperation in following an appropriately defined postoperative management program, and conducting postoperative follow-up examinations.

**Preoperative:**

- As part of the preoperative examination, the surgeon must check that no factors, especially biological and biomechanical, will affect the correct performance of the implant during the operation and postoperative period.
- Only patients that meet the criteria described in the indications For Use should be selected.
- Patient conditions and/or pre-dispositions such as those addressed in the aforementioned contraindications should not be selected.
- Care should be used in the handling and storage of the Minuteman components. Implants and instruments should be protected during storage especially from corrosive environments.
- The type of construct to be assembled for the case should be determined prior to beginning the surgery. Based on the fatigue testing results, the surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may have an impact on the performance of the Minuteman. An adequate inventory of implant sizes should be available at the time of surgery, including an excess of implant sizes expected to be used and implant sizes larger and smaller than those expected to fit in the patient.
- Since mechanical parts are involved, the surgeon should be familiar with the various components before using the equipment, and should personally assemble each device to be used to verify that all parts and necessary instruments are present before the surgery begins. Damaged or defective instruments should not be used. Contact the manufacturer for repair or replacement instructions.
- All components should be cleaned and sterilized before use. Additional sterile components should be available in case of an unexpected need.

**Intraoperative:**

- At all times, extreme caution should be used around the spinal cord and nerve roots. Damage to nerves may result in loss of neurological functions.
- Breakage, slippage or misuse of instruments or implant components may cause injury to the patient or operative personnel.
- The correct selection of the type and size of implant appropriate to the patient, and positioning of the implant, are extremely important to allow maximum chances for a successful surgical result.
- Bone grafts, if used, must be placed in the area to be fused and be in contact with viable bone.

**Postoperative:**

- The physician's postoperative directions and warnings to the corresponding patient are important to allow maximum chances for a successful surgical result. Detailed instructions on the use and limitations of the device should be given to the patient. If excessive or early weight-bearing, or excessive muscular activity is discouraged during the early postoperative rehabilitation period, the patient must be warned that bending, loosening or breakage of the components are complications which can occur as a result of this activity. The risks of bending, loosening or breakage of the components are complications which can occur as a result of this activity, or if the patient is debilitated, demented, or using weight supporting devices. The patient should be warned to avoid falls or sudden jolts to lessen the possibility for bending, loosening or breakage of the internal fixation device.
- To allow maximum chances for a successful surgical result, the patient or device should not be exposed to mechanical vibrations that may loosen the device construct. The patient should be warned of this possibility, and instructed to limit physical activities, especially lifting, twisting and any type of sport participation. The patient should be advised not to smoke, utilize nicotine products or consume alcohol or nonsteroidal anti-inflammatory drugs such as aspirin during the bone graft healing process.
- The patient should be advised of their inability to bend at the point of spinal fusion and taught to compensate for this permanent physical restriction in body motion.
- If a nonunion develops or if the components loosen, bend, and/or break, the device(s) should be revised and/or removed immediately before serious injury occurs. Failure to immobilize a delayed nonunion of bone will result in excessive and repeated stresses on the implant. By the mechanism of fatigue these stresses can cause eventual bending, loosening, or breakage of the device(s).
- The patient should be fully instructed in the appropriate postoperative care. The patient's ability and willingness to follow, as well as comprehension of the importance of following, instructions are one of the most important aspects of successful postoperative healing.
- Explanted surgical implants must never be reused.
- As a precaution, before patients with implants receive any subsequent surgery (such as dental procedures), prophylactic antibiotics should be considered, especially for patients with increased risk for infection.

**Sterility:**

Gamma irradiation is indicated by the "Sterile-R" symbol on the Minuteman implant labeling. These devices remain sterile as long as the package integrity has not been violated.

**Packaging:**

Packages for each of the Minuteman components should be intact upon receipt. Inspect each package prior to use and do not use the component if any seal or cavity is damaged or breached or if the expiration date has been exceeded. Once opened, the component must be used, discarded, or returned to manufacturer. If a loaner or consignment instrument system is used, all sets should be carefully checked for completeness to ensure there is no damage prior to use. Damaged packages or products should never be used and should be returned to the manufacturer.

**Handling and storage:**

The Minuteman implants and instruments must be stored with care. Before use, inspect all instrumentation for proper function, possible damage, wear or non-function. Damaged or defective instruments should not be used. Note: that at some point in time, instruments wear out and should be replaced. Contact the manufacturer for repair or replacement instructions.

**Recommendations for Care and Cleaning of Minuteman Surgical Instrument System:**

Both physical and chemical (detergent) processes are necessary to minimize the bioburden on all soiled items. Chemical (detergent) cleaners alone cannot remove all soil and debris, therefore a careful manual cleaning of each item is essential for maximum decontamination. Spinal Simplicity recommends the use of a mild enzymatic detergent with a near-neutral pH prepared following the manufacturer's instructions for preparation and use. Saline solution should NOT be used, as saline has a corrosive effect on stainless steel. Remove heavy or large debris using single-use, non-shedding wipes soaked in cleaning solution. Immerse instruments in prepared bath. Using a soft bristle brush (do not use steel brushes) brush all surfaces of the instruments while they are submerged in bath, ensuring that all visible soil is removed. Whenever applicable: use a pipe cleaner and syringe to clean all cannulae, lumens, crevices, grooves and hard to reach areas. Repeatedly operate/bend/articulate all moveable components while cleaning. Allow instruments to soak in detergent prepared bath for manufacturer's recommended soaking time. Final cleaning and disinfecting should be performed via sonication or an automated washer/disinfectant cycle. See below recommendations for each:

**Sonication:**

1. Remove the instruments from bath and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumen, holes, and difficult-to-reach areas.
2. Place prepared cleaning agents in a sonication unit. Completely submerge instruments in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
3. Rinse instruments in purified water for at least 3 minutes.
4. Visually inspect to determine if all visible soil has been removed from the surfaces, lumen, cannulae, crevices, serrations, threading, etc. If visible soil remains, repeat the cleaning/disinfecting procedure.
5. Dry the instruments with single-use, non-shedding absorbent wipes and/or medical compressed air (e.g. interiors and cannulae). Be sure to completely dry the instruments immediately after rinse to inhibit corrosion.

**Automated Washer:**

1. Remove instruments from bath and rinse in purified water for a minimum of 1 minute. Thoroughly and aggressively flush lumen, holes, and difficult-to-reach areas.
2. Place instruments in a suitable washer/disinfectant basket and process through a standard washer/disinfectant cleaning cycle. The following minimum parameters are essential for thorough cleaning and disinfecting.
  - a. 2 minute prewash with cold tap water
  - b. 20 second enzyme spray with hot tap water
  - c. 1 minute enzyme soak

- d. 15 second cold tap water rinse
- e. 2 minute detergent wash with hot tap water (64-66°C/146-150°F)
- f. 15 second hot tap water rinse
- g. 2 minute thermal rinse (80-93°C/176-200°F)
- h. 10 second purified water rinse (64-66°C/146-150°F)
- i. 7 to 30 minute hot air dry (116°C/240°F)

NOTE: The washer/disinfectant manufacturer's instructions should be strictly adhered to.  
3. Visually inspect to determine if all visible soil has been removed from the surfaces, lumen, cannulae, crevices, serrations, threading, etc. If visible soil remains, repeat the cleaning/disinfecting procedure.

Caution: certain cleaning solutions such as those containing formalin, glutaraldehyde, mercury, active chlorine, chloride, bromine, bromine, iodine, and/or alkaline cleaners may damage some instruments. Such cleaning solutions should not be used.

#### Recommendations for Sterilization of Minuteman Surgical Instrument System:

For typical steam autoclave cycles, the following are recommended times and temperatures:

1. **Prevacuum Sterilizer:**  
Wrapped cases, trays and instruments, or cases, trays and instruments should be exposed to 135° C (275° F) for at least 3 minutes. Dry for 30 minutes.

#### Product Complaints:

Communicate suspected deficiencies in product quality, identity, durability, reliability, safety, effectiveness and/or performance directly to Spinal Simplicity, LLC. When filing a complaint provide component name(s), part number(s), lot number(s), your name and address, the nature of the complaint and patient case number. Sterilize and return all component(s) to Spinal Simplicity, LLC. Notify Spinal Simplicity, LLC immediately of any complaints and any incidents resulting in patient death or serious injury.

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